Home > ... > FP7 >

Macromolecular Dynamics in an Optimized Microfluidic Cross-Slot Geometry: Application to Biofluid Analogues and Prosthetic Fluid Formulation





Macromolecular Dynamics in an Optimized Microfluidic Cross-Slot Geometry: Application to Biofluid Analogues and Prosthetic Fluid Formulation

Results in Brief

Device to detect changes in joint fluid

Altered rheology of bodily fluids is associated with numerous pathological conditions. A microfluidic device to measure viscosity of very small samples could pave the way to rapid diagnosis, monitoring and even therapy for osteoarthritis.





[©] Thinkstock

Some materials are sort of stretchy and stringy due to the presence of long and flexible molecules that resist deformation under an applied longitudinal stress. Melted polymers and melted cheese are such materials. So are important biological materials such as saliva, mucus, and the synovial fluid that lubricates and protects cartilage surfaces in mammalian joints.

In arthritis patients, the viscoelasticity of synovial fluid is reduced. A microfluidic platform (a miniature flow channel) developed within the EU-funded project MICRO-BIORHEOLOGY can now characterise synovial fluid undergoing stretching deformations. It could be applied to other fluids as well.


Scientists focused on hyaluronic acid, a polysaccharide found in large quantities in synovial fluid. Hyaluronic acid injections are one treatment for osteoarthritis. Improved understanding will lead to both better diagnosis based on fluid samples and better treatments.

The microfluidic extensional flow device, termed optimised-shape cross-slot extensional rheometer (OSCER), exploits mutually bisecting channels (cross-slots) in order to stretch microscopic fluid samples. The team first confirmed that it accurately reproduces computational fluid dynamics simulations of the flow of simple fluids like water and 'model' viscoelastic fluids of polymers. Using polarised light microscopy on polymer solutions, scientists showed that changes in the intensity of transmitted light are proportional to fluid stress. This makes these changes useful in measuring resistance to extensional deformation (i.e. the extensional viscosity).


Investigators then measured the extensional viscosity of model synovial fluids based on hyaluronic acid solutions. For each fluid tested, the OSCER device measured extensional viscosity up to 50 times higher than the shear viscosity measured by a conventional rheometer. This shows that the synovial fluid 'hardens' in response to strain. In practice, this would mean that it thickens as it is compressed between the joint surfaces, providing shock absorption and preventing damage, particularly from sudden high-load impacts.

Further studies demonstrated that refractive index differences can be used to assess hyaluronic acid molecular properties. This supports the possibility of characterising hyaluronic acid in real synovial fluid to evaluate degradation. The methods could provide the basis for a diagnostic toolbox that could also follow the progression of the disease with biopsy-sized samples of synovial fluid.

Keywords



Discover other articles in the same domain of application



Grant agreement ID: 298220		Specific programme "People" implementing the
		Seventh Framework Programme of the European
Project closed		Community for research, technological
		development and demonstration activities (2007 to
		2013)
Start date	End date	2013)
2 July 2012	1 July 2014	
		Total cost
		€ 208 672,40
		EU contribution
		€ 208 672,40
		Coordinated by
		UNIVERSIDADE DO PORTO
		Portugal

Last update: 13 April 2015

Permalink: <u>https://cordis.europa.eu/article/id/159642-device-to-detect-changes-in-joint-fluid</u>

European Union, 2025