Home > ... > H2020 >

Proof of concept of prognostic test for critical Septic Patients based in Early Detection of NLRP3 Inflammasome impairment activation.



Proof of concept of prognostic test for critical Septic Patients based in Early Detection of NLRP3 Inflammasome impairment activation.

Results in Brief

Point-of-care device quickly detects early signs of sepsis

With sepsis killing more than 5 million people each year globally, SPEDI-TEST has developed a mobile device to quickly detect which patients are likely to develop complications – leading to early life-saving interventions.





© PhotoByToR/Shutterstock.com

Sepsis occurs when the body's immune system overreacts to infection or injury and ends up attacking its own organs and tissues. In 2017 there were estimated to be almost 50 million cases globally

Despite the availability of prediction and diagnostic models, clinicians are still unable to assess the gravity of sepsis in a patient or the likely outcome.

The SPEDI-TEST project, supported by the <u>European Research Council</u> As developed a prognostic method, which detects sepsis within hours.

"Our device analyses specific blood cell culture conditions to spot when the

inflammatory response is impaired, then assesses inflammation at the molecular level, as a biomarker for sepsis-related late complications," says Pablo Pelegrin, project coordinator.

The team has launched a <u>spin-off company</u>, to convert their in vitro diagnostic (IVD) device into an automated compact lab-on-a-chip device.

Identifying the blood biomarkers

In sepsis, the inflammatory response is often triggered by proteins called receptors, which sit on the surface of immune cells. Once activated, these secrete proinflammatory cell signalling proteins called cytokines that travel around the body activating immune cells to fight infection.

SPEDI-TEST developed two protocols for identifying blood biomarkers for sepsis.

For the first, different patterns of specific proteins secreted by the cells of patients with sepsis are detected in purified blood samples using <u>ELISA</u>. From this observation, the team then used flow cytometry to study each patient's immune system, enabling predictions about how sepsis is likely to develop.

In the second protocol, complete blood samples are used (instead of purified blood cells) with specific cell culture conditions. While this will make future automation easier, whole blood samples are less sensitive than purified cells. The team again used both ELISA and flow cytometry, finding flow cytometry more sensitive, with results still undergoing analysis.

"Rather than looking for a molecule present in the blood, we 'ask' the immune cells how they are doing, in patients with long-term complications; they effectively answer, 'we are exhausted', which becomes our biomarker," explains Pelegrin.

This sort of blood work usually requires cumbersome equipment and complicated techniques and specialists.

"Our point-of-care IVD device is small enough to be taken to the patient's bedside and can be used easily by nurses or clinical staff, without special training. Its fast turnaround of results, in a few hours, speeds up potentially life-saving interventions," adds Pelegrin.

Validation outcomes

Validation with patients found ELISA to have a sensitivity rate (correctly identifying patients who will develop future sepsis complications) of 90.9 %, and a specificity

(correctly identifying those without) of 100 %. Whereas the flow cytometry achieved a sensitivity of 81.8 % and specificity of 91.7 %. Both detected patients with long-term sepsis-related complications within 24 hours of admission to intensive care units.

"This should be compared to the current clinical measures, known as Δ SOFA, used to determine likely sepsis evolution calculated five days after ICU admission of 55 patients, which has sensitivity and specificity rates of 87.5 %," adds Pelegrin.

Better identification of risk, more effective treatment

SPEDI-TEST identifies early patients with an increased chance of developing lifethreatening sepsis complications, resulting in earlier prognosis and so more effective treatment and care. It also benefits vulnerable populations, for example those with comorbidities such as myocardial infarctions or pneumonia.

"With the cost of care per septic patient estimated at EUR 42 000 annually, this benefits not only patients, but also hospitals and the wider healthcare economy through more targeted resources," concludes Pelegrin.

Keywords

SPEDI-TEST, sepsis, immune system, diagnostic, biomarker, inflammatory, blood, ELISA, cytometry, cytokines, molecule, protein, IVD

Discover other articles in the same domain of application



Are there really bacteria in the womb?






AX-1 the farmer robot conquers weeds selectively and sustainably





When you do not speak or sign the same language





Investigating blood flow in severe COVID-19 patients



Project Information

SPEDI-TEST

Grant agreement ID: 899636

DOI 10.3030/899636

Project closed

EC signature date 16 June 2020

Funded under EXCELLENT SCIENCE - European Research Council (ERC)

Total cost No data

EU contribution € 150 000,00

Coordinated by

Start date 1 October 2020 End date 31 March 2022 FUNDACION PARA LA FORMACION E INVESTIGACION SANITARIAS DE LA REGION DE MURCIA Spain

Related articles



Last update: 21 October 2022

Permalink: <u>https://cordis.europa.eu/article/id/442330-point-of-care-device-quickly-detects-early-signs-of-sepsis</u>

European Union, 2025